

Implant Hybrid Bridge



Precision milled
titanium bar

Close fit between
the two parts

Esthetic finishes for
both teeth and gum



Choices of two materials
for the Superstructure

1. High Translucency Zirconia

2. High Strength Acrylic

Features

- 1. Titanium bar provides**
 - i. Strong connections to implants
 - ii. Full support for super structure
- 2. Thick super structures more resistant to breakage**
- 3. Long distal extension possible**
- 4. Full digital workflow with Photogrammetry device**

Digital Workflow



1. DATA ACQUISITION



2. PHOTOGRAMMETRY



**Tupel 3D Implant Scanner*



3. CAD-DESIGN



4. CAM-PRODUCTION



5. ESTHETIC FINISH



***Contact Us
For More Details**



Full Digital Workflow



**Strong Connections
to Implants**



**Support from
Titanium Bar**



**Choices of
Superstructure**



Affordable Price